## medilog®AR

### Efficient and reliable state-of-the-art Holter



# medilog®AR

- Robust, shock- and splashproof as well as easy to clean
- Flexible dual-battery concept for more than 14 days recording duration
- Zero seconds atrial fibrillation/atrial flutter detection based on true P-wave analysis

#### IDEALLY SUITED FOR CLINICAL ROUTINE

The shock- and splash-proof design makes the medilog®AR Holter highly suitable for the daily clinical routine. The Holter and its carrying case are robust and easy to clean.

#### **TWO WEEKS OF POWER**

Thanks to the dual-battery concept, the patients can be screened for more than 14 days without having to see the doctor to change the batteries.



Eco-friendly dual-battery concept (illustration)

#### **ECO-FRIENDLY BATTERY:**

By adjusting the signal resolution, the internal battery can run for 24/72 hours or more than 7 days. The high-quality internal battery is recharged via USB.

#### **AAA SINGLE USE BATTERY:**

For fastest turnaround times in high throughput clinics, use easily available standard AAA batteries.

#### **MORE THAN 14 DAYS:**

For maximum recording duration, combine a fully charged internal battery with an AAA battery to get more than 14 days with automated AF detection.

#### VARIOUS CARRYING OPTIONS

medilog<sup>®</sup>AR is so flat that it can easily be carried in a trouser pocket. With the carrying case, it can be worn on a belt or with a neck belt directly on the body.

#### MADE TO SUPPRESS ARTIFACTS

The medilog®AR Holter recorder features a superior resting ECG grade 32,000Hz sampling rate on 3 channels. This is used for true P-wave analysis, sophisticated artifact suppression and motion detection, enabling a speedy evaluation.

#### ZERO SECONDS ATRIAL FIBRILLA-TION/AFL DETECTION:

The medilog®AR scans for atrial fibrillation and atrial flutter by direct-ly looking for P-waves, just as a cardiologist would.

#### display for easy operation

robust carrying case with various carrying options microphone to record patient ID

buttons for patient events

runs on an internal rechargeable battery, or on normal AAA batteries, if desired

plug and play USB data exchange

removeable patient cable with 5 wires, alternatively with 7 wires for even more reliability

medilogA

#### **RESPIRATION EPISODES:**

Thanks to ECG-derived respiration recording, the medilog®AR is able to screen for potential respiratory episodes during sleep. The optional SpO<sub>2</sub> sensor, connected via Bluetooth, allows for additional respiratory information.



 $\mathsf{SpO}_2$  sensor, attached for the night

#### MUCH MORE THAN ARRHYTHMIA SCREENING

Locating Atrial Fibrillation in zero seconds, analysis of respiratory episodes, and high-resolution HRV are just some of the exceptional features of the medilog®AR. The analysis is carried out using the renowned medilog®DARWIN2 software.



Americas SCHILLER Americas Inc. Doral, Florida 33172

North America: Phone +1 786 845 06 20 Fax +1786 845 06 02 sales@schilleramericas.com www.schilleramericas.com

Latin America & Caribbean Phone +1 305 591 11 21 Fax +1 786 845 06 02 sales@schilleramericas.com www.schilleramericas.com



Asia SCHILLER Asia-Pacific / Malaysia 52200 Kuala Lumpur, Malaysia Phone +603 6272 3033 Fax +603 6272 2030 sales@schiller.com.my www.schiller-asia.com



Austria SCHILLER Handelsgesellschaft m.b.H. A-4040 Linz Phone +43 732 70 90 Fax +43 732 757 000 sales@schiller.at www.schiller.at







10000 Zagreb Phone +385 1 309 66 59 Fax +385 1 309 66 60 info@schillerzg.hr www.schiller.hr



France SCHILLER Médical F-67160 Wissembourg Phone +33 3 88 63 36 00 +33 3 88 63 36 49 Fax info@schiller.fr www.schiller-medical.com



France (distribution France) SCHILLER France S.A.S. F-77608 Bussy St Georges Phone +33164 66 50 00 Fax +33164665010

contact@schillerfrance.fr

www.schiller-france.com



Germany SCHILLER Medizintechnik GmbH D-85622 Feldkirchen b. München Phone +49 89 62 99 81 0 Fax +49 89 609 50 90 info@schillermed.de www.schillermed.de



SCHILLER Healthcare India Pvt. Ltd. Mumbai – 400 059, India

Phone +91 22 6152 3333 / 2920 9141 Fax +91 22 2920 9142 sales@schillerindia.com www.schillerindia.com



schiller@schiller.pl www.schiller.pl Russia & C.I.S.

AO SCHILLER.RU 125040 Moscow, Russia Phone +7 (495) 970 11 33 Fax +7 (495) 956 29 10 mail@schiller.ru www.schiller.ru



SCHILLER d.o.o. 11010 Beograd Phone +381 11 39 79 508 Fax +381 11 39 79 518 info@schiller.rs www.schiller.rs



Slovenia SCHILLER d.o.o. 2310 Slovenska Bistrica Phone +386 2 843 00 56 Fax +386 2 843 00 57 info@schiller.si www.schiller.si



SCHILLER ESPAÑA, S.A. ES-28232 Las Rozas/Madrid Phone +34 91 713 01 76 Fax +34 91 355 79 33 schiller@schiller.es www.schiller.es



SCHILLER-Reomed AG CH-8912 Obfelden Phone +41 44 744 30 00 Fax +41 44 740 37 10 sales@schiller-reomed.ch www.schiller-reomed.ch



Turkey SCHILLER TÜRKIYE Okmeydani-Sisli – Istanbul Phone +90 212 210 8681 (pbx) Fax +90 212 210 8684 info@schiller.com.tr www.schiller-turkiye.com



SCHILLER UK Bellshill, ML4 3PR Phone +44 1698 744 505 Fax +44 1698 744 474 sales@schilleruk.com www.schilleruk.com

Device availability in your market is subject to regulatory approval.



**SCHILLER** The Art of Diagnostics